


Piezo Guide

For the guide mechanism, Sigma fine stages adopted a guide system that utilizes elastic deformation of metals and a mechanism to increase deformation of piezo elements. These originally designed stages achieved readable resolution of 10nm during closed loop control, ideal for uses that require high-speed high-precision positioning.


Structure of Sigma Fine Stage

Displacement Magnification Method

Piezo actuator and displacement magnification mechanism offer a large operating range.



SFS-H Internal Structural Drawing



Digital Sensor

Closed loop control is possible with a digital sensor that does not require any high precision analog amplifier or AD conversion circuit.

Operating Environment of Sigma Fine Stage


Use fine stages within the following operating environment temperature range. Contact us separately if you desire to use the stages outside the operating environment temperature range.

***Operating environment**
 Temperature: 10°C – 30°C
 Humidity: 20% – 60%
 (without condensation)

***Recommended environment**
 Temperature: 20°C ±1°C
 Humidity: 40% or lower

Since durability of piezo elements used in the SFS/SFS-H stage series will deteriorate in high humidity environments, use them in the above environments.

Sigma Fine System Chart



- Application Systems
- Optics & Optical Coatings
- Holders
- Bases
- Manual Stages
- Actuators

Motoeized Stages

- Light Sources
- Index

- Guide
- Controllers/Drivers
- Softwares
- Stepping Motor
- AC Servo Motor
- Cables


Piezo

- X Translation
- Theta Rotation
- Goniometer
- Vacuum
- Options


- 40mm
- 60mm
- 80mm
- 85mm
- 100mm
- 120mm
- Others

Accuracy Measurement Example: Sigma Fine Stage System SFS-H (Linear Stage)

High Speed Convergence Data




0⇔75µm pulse rate (16Hz)
Closed loop control
by SFS-H60X at 150g load




0⇔1µm step convergence data (15msec)
Closed loop control
by SFS-H60X at 150g load

Travel

The following graph shows the hysteresis curve unique to piezo actuators during open loop control travel.




Linearity



0⇔80µm linearity within 0.3%
Closed loop control
by SFS-H60X at 150g load


Positional Repeatability



0⇔80µm positional repeatability 50nm or lower
by SFS-H60X at 150g load


Characteristics of Fine Feed

Characteristics when feed amount is small in closed loop control. Hysteresis disappears in open loop control.




Straightness

Deviation from the straight line in the direction of travel.




Pitch

Tilt around the axis in the horizontal plane perpendicular to the direction of travel.



Characteristics of Yaw


Rotation around the axis in the vertical plane perpendicular to the direction of travel.



Follow-up example with Respect to Analog Input: SFS-H (Linear Stage) *Controlled by FINE-01γ


High Speed Convergence Data

Input waveform: Navy/Output waveform: Light blue




SFS-H40X

Staircase wave input/output waveforms
(Step 50nm 20-step staircase wave)




SFS-H40X

Saw-tooth wave input/output waveforms
(Uniform motion 35µm/1Hz)



SFS-H40X

Rectangular wave input/output (30Hz)



SFS-H40X

Sine wave input/output (30Hz)

Sigma Fine (Piezo) Stages (high stiffness type) XY Piezo Stages Aperture Type

SFS-H
SFS-120XY(WA)



These piezo stages achieved fine movement of theoretical resolution of 1nm, offering high precision, high stiffness and high speed with adoption of digital sensors.

SFS-H




- These piezo stages offer high precision and high resolution positioning because full closed loop control is possible with digital sensors.
- With piezo elements used as actuators, travel between 90µm – 100µm in open loop control, and fine movement of theoretical resolution of 1nm are possible.
- Readable resolution of 10nm is achieved during full closed loop control with adoption of a digital sensor (micro-displacement sensor based on the frequency digital conversion method) used as position detection sensors.
- Sigma Koki controller FINE series offers high speed. [Reference](#) G117 Regarding pulse rate, SFS-H (linear) achieved 10 – 15Hz, or 25Hz when the current FINE-01γ analog input is used.




Outline Drawing


SFS-H40X(CL) Hexagon socket head cap screw M3×10...2 screws




SFS-H40Z(CL) Hexagon socket head cap screw M3×10...2 screws



SFS-H60X(CL) Hexagon socket head cap screw M4×10...4 screws



SFS-H60XY(CL) Hexagon socket head cap screw M4×10...4 screws



Specifications

| Part Number | SFS-H40X(CL) | SFS-H40Z(CL) | SFS-H60X(CL) | SFS-H60XY(CL) |
|--|---------------------|---------------------|--|--|
| Travel | 90µm±15% | 100µm±15% | 100µm±15% | 100µm±15% |
| Table Size [mm] | 40×40 | 40×40 | 60×60 | 60×60 |
| Actuator | Piezo actuator | Piezo actuator | Piezo actuator | Piezo actuator |
| Weight [kg] | 0.28 | 0.28 | 0.4 | 0.43 |
| Theoretical Resolution (open-loop) [nm] | 1 | 1 | 1 | 1 |
| Resolution (closed-loop) [nm] | 10 | 10 | 10 | 10 |
| Linearity [%] | 0.3 or lower | 0.3 or lower | 0.3 or lower | 0.3 or lower |
| Perpendicularity (Horizontal Direction) [µm] | 1 | 1 | 1 | 1 |
| Positional Repeatability [µm] | 0.1 or lower | 0.1 or lower | 0.1 or lower | 0.1 or lower |
| Load Capacity [N] | 9.8 (1.0kgf) | 6.7 (0.7kgf) | 19.6 (2.0kgf) | 14.7 (1.5kgf) |
| Micro-displacement Sensor | Digital Sensor | Digital Sensor | Digital Sensor | Digital Sensor |
| Compatible Cable | Attached cable (2m) | Attached cable (2m) | FINE-CA-3: For piezo DS1-CA-3: For digital sensor | FINE-CA-3: For piezo DS1-CA-3: For digital sensor |

Application Systems
Optics & Optical Coatings
Holders
Bases
Manual Stages
Actuators

MotORIZED Stages

Light Sources

Index

Guide

Controllers/Drivers

Softwares

Stepping Motor

AC Servo Motor

Cables

Piezo

X Translation

Theta Rotation

Goniometer

Vacuum

Options

40mm
60mm

80mm
85mm


100mm
120mm

Others




Outline Drawing

SFS-H60Z(CL) Hexagon socket head cap screw M4x10...4 screws




SFS-H60XYZ(CL) Hexagon socket head cap screw M4x10...4 screws




SFS-120XY(WA)

- A high precision positioning XY axis stage with 50x50mm aperture fitted with a piezo element as an actuator and digital sensor for position feedback, suitable for fine movement of samples under a microscope.
- Can be driven with the Sigma fine stage controller FINE-503. Since the controller supports RS232C, GP-IB (FINE-503 only) and USB interfaces, position control can be performed easily from a PC using the software for positioning & measurement SGEMCSE, SGTERME and SGSFSXE.
- Ideal for applications such as incorporation into various types of microscopes, precision measurement, semiconductor test equipment, high-precision mask alignment, scanning interferometer, image processor, and biotechnology systems.



Outline Drawing

SFS-120XY(WA) Hexagon socket head cap screw M4x8...4 screws



Specifications

| Part Number | SFS-H60Z(CL) | SFS-H60XYZ(CL) | SFS-120XY(WA) |
|--|--|--|--|
| Travel | 100μm±15% | 100μm±15% | 100μm±10% |
| Table Size [mm] | 60x60 | 60x60 | 120x120 |
| Actuator | Piezo actuator | Piezo actuator | Piezo actuator |
| Weight [kg] | 0.33 | 0.63 | 1.2 |
| Theoretical Resolution (open-loop) [nm] | 1 | 1 | 1 |
| Resolution (closed-loop) [nm] | 10 | 10 | 10 |
| Linearity [%] | 0.3 or lower | 0.5 or lower | — |
| Perpendicularity (Horizontal Direction) [μm] | 1 | 1 | 1 or lower |
| Positional Repeatability [μm] | 0.1 or lower | 0.15 or lower | 0.1 or lower |
| Load Capacity [N] | 9.8 (1.0kgf) | 9.8 (1.0kgf) | 19.6 (2.0kgf) |
| Micro-displacement Sensor | Digital sensor | Digital sensor | Digital sensor |
| Compatible Cable | FINE-CA-3: For piezo DS1-CA-3: For digital sensor | FINE-CA-3: For piezo DS1-CA-3: For digital sensor | FINE-CA-3: For piezo DS1-CA-3: For digital sensor |

Application Systems

Optics & Optical Coatings

Holders

Bases

Manual Stages

Actuators

MotORIZED Stages

Light Sources

Index

Guide

Controllers/Drivers

Softwares

Stepping Motor

AC Servo Motor

Cables

Piezo

X Translation

Theta Rotation

Goniometer

Vacuum

Options

□40mm

□60mm

□80mm

□85mm

□100mm

□120mm

Others

Piezo Actuator for Objective Lens

SFS-OBL (Upright)/SFAI-OBL (Inverted)



Catalog Code **W9056**

Objective lens actuator for inverted microscope employing a piezo element as actuator and digital sensor for feedback.




- Compact, and enabling high-speed high-resolution positioning.
- Travel is 100µm at open loop.
- Two types of erected model and inverted model are available for incorporation into various types of microscopes.
- In the case of the Sigma fine stage series, these actuators can be driven with the controller (FINE-01γ/503(CL)). Since RS232C, GP-IB (FINE-503 only) and USB interfaces are supported, position control can be performed easily from a PC using the software for positioning & measurement SGEMCSE, SGTERME and SGSFSXE.

Guide

- ▶ Adapters compatible with screw sizes of other manufacturers' objective lenses are also available.
- ▶ SFS-OBL-2 uses a metal enclosure type piezo actuator to improve environment resistance such as humidity compared to SFS-OBL-1.


Outline Drawing



Specifications

| Part Number | CE SFS-OBL-1 | SFS-OBL-2 | SFAI-OBL-1R |
|---|-------------------------|-------------------------|---------------------------------------|
| Travel | 100µm±15% | 100µm±15% | 100µm±15% |
| Objective Lens Diameter [mm] | Diameter φ39 or less | Diameter φ39 or less | Diameter φ39 or less |
| Dimensions [mm] | (W)75.5 × (H)45 × (D)40 | (W)75.5 × (H)55 × (D)40 | (W)60.8 × (H)30 × (D)40 |
| Actuator | Piezo element | Piezo element | Piezo element |
| Weight [kg] | 0.15 | 0.24 | 0.15 |
| Theoretical Resolution (open-loop) [nm] | 1 | 1 | about 0.8 |
| Resolution (closed-loop) [nm] | 10 | 10 | 10 |
| Straightness (Xy Xz Yx Yz) [µm] | 1 or lower | 1 or lower | 0.2 or lower |
| Positional Repeatability [µm] | 0.1 or lower | 0.1 or lower | 0.1 or lower |
| Load Capacity [N] | — | — | 4.9 (0.5kgf) |
| Micro-displacement Sensor | Digital sensor | Digital sensor | Digital sensor |
| Compatible Adapter | OBL-ADP-** | OBL-ADP-** | OBL-ADP3-** |
| Accessories | Cable (2m) | Cable (2m) | Cable (2m), four special lift spacers |

Objective Lens Adapters



SFS-OBL Compatible Adapters

| Part Number | Mounting Screw Size [mm] | A [mm] | B [mm] |
|-----------------|----------------------------|----------------------------|--------|
| OBL-ADP-M20.32A | Microscope side M20.32 | M20.32 P=0.706 (W0.8×1/36) | 15 |
| OBL-ADP-M20.32B | Objective lens side M20.32 | M20.32 P=0.706 (W0.8×1/36) | — |
| OBL-ADP-M25.0A | Microscope side M25.0 | M25.0 P=0.75 | 20 |
| OBL-ADP-M25.0B | Objective lens side M25.0 | M25.0 P=0.75 | — |
| OBL-ADP-M26.0A | Microscope side M26.0 | M26.0 P=0.706 (W26.0×1/36) | 21 |
| OBL-ADP-M26.0B | Objective lens side M26.0 | M26.0 P=0.706 (W26.0×1/36) | — |

SFAI-OBL Compatible Adapters

| Part Number | Mounting Screw Size [mm] | A [mm] | B [mm] |
|------------------|----------------------------|----------------------------|--------|
| OBL-ADP3-M20.32A | Microscope side M20.32 | M20.32 P=0.706 (W0.8×1/36) | 15 |
| OBL-ADP3-M20.32B | Objective lens side M20.32 | M20.32 P=0.706 (W0.8×1/36) | — |
| OBL-ADP3-M25.0A | Microscope side M25.0 | M25.0 P=0.75 | 20 |
| OBL-ADP3-M25.0B | Objective lens side M25.0 | M25.0 P=0.75 | — |
| OBL-ADP3-M26.0A | Microscope side M26.0 | M26.0 P=0.706 (W26.0×1/36) | 21 |
| OBL-ADP3-M26.0B | Objective lens side M26.0 | M26.0 P=0.706 (W26.0×1/36) | — |

- Application Systems
- Optics & Optical Coatings
- Holders
- Bases
- Manual Stages
- Actuators

MotORIZED Stages

Light Sources

Index

Guide

Controllers/Drivers

Softwares

Stepping Motor

AC Servo Motor

Cables

Piezo

X Translation

Theta Rotation

Goniometer

Vacuum

Options

- 40mm
- 60mm
- 80mm
- 85mm
- 100mm
- 120mm

Others

Controllers with built-in piezo drivers for single axis / 3 axes.



FINE-01y



FINE-503(CL)

- These controllers are fitted with digital sensor input for each axis, enabling closed loop control by correcting hysteresis curve unique to piezo.
- Being connected to a PC via RS232C, GP-IB, or USB interface, the FINE-503(CL) allows control of a fine stage by simple commands sent from a PC.
- In addition to PC control via RS232C or USB interface, the FINE-01y enables high-speed control synchronized with analog signal input.

| Part Name | Part Number |
|--|----------------------|
| 1 axis SFS Controller with Analog Input Function | FINE-01y |
| 3 axes SFS Controller | FINE-503(CL) |
| Control Pad | CJ-200A |
| FINE Cable | FINE-CA-3 |
| DS Cable | DS1-CA-3 |
| BNC-BNC Cable | SKBNC-BNC-3.0 |

Primary Functions

| Part Number | FINE-01y | FINE-503(CL) |
|------------------------|----------------|--------------|
| Controller Function | | ○ |
| Number of Control Axes | 1 | 3 |
| Stored Program Control | | ○ |
| Feedback Control | Digital sensor | |

General Specifications

| | |
|-------------------------------|-------------------------|
| Power Voltage | AC100V ±10% 50/60Hz |
| Power Consumption | 50VA |
| Operating Temperature | 10 – 30°C |
| Storage Temperature | -20 – 60°C |
| Ambient Humidity | 20 – 80%RH |
| External Dimensions (WxHxDmm) | 225x118x250/270x118x297 |
| Weight (kg) | 3.5 5.3 |

Interface

| | | |
|--------------|---|---|
| GP-IB | — | ○ |
| RS232C | | ○ |
| USB | | ○ |
| Analog input | ○ | — |

Optional

| | | |
|---------------|---|---|
| CJ-200A | | ○ |
| SKBNC-BNC-3.0 | ○ | — |

Performance Specifications

| | |
|-----------------------------|------------|
| Coordinate Indication Range | ±999,999nm |
| Max. Travel to Set | ±999,999nm |


Control Command

| | |
|--------------------------------|---|
| Machine Origin Return | ○ |
| Theoretical Origin Setting | ○ |
| Relative Position Drive | ○ |
| Absolute Position Drive | ○ |
| Jog Operation | ○ |
| Position Appointment | — |
| Circular Interpolation Control | — |
| Linear Interpolation Control | — |
| Drive | ○ |
| Deceleration Stop | — |
| Emergency Stop | — |
| Speed Setting | ○ |
| Motor Free/Hold | — |
| Port Input | — |
| Port Output | — |


SFS Software

Free Software | SGCommander (for RS232C) Windows® Version

Free software is available to operate your controller easily from a PC. Each axis of a connected motorized stage can be moved using buttons on the screen. The software can be downloaded from our website.




Simple operations are possible such as travel by specifying an axis, homing or jog operation.



Controllers such as SHOT-30*/702 and FINE-**, which have a built-in program function, allow editing of programs from a PC. Since data can be downloaded/uploaded from/to Excel sheets, it is easy to edit programs. In addition, upload of memory switch or download mode is available.

Free Application LabVIEW (for v.5.1/v.6i/v.7.1/v.8.6/v.2010) RS232C/GP-IB

LabVIEW application is available for LabVIEW users.



Other: 30 Day Trial Version (SGEMCSE/SGTERME)

SGEMCSE....SGEMCSE is software for collecting data or measuring using automatic positioning equipment, measuring instrument or controller, and is offered 30 days for free.






SGTERME....It allows command input using Excel for easy program making. SGTERME is ideal software to link with various devices, and is offered 30 days for free.



With Industrial Line positioners, attocube has genuinely combined highest precision piezodrive technology with extremely rugged, yet cost effective design.

- All ECS positioners of the Industrial Line are dedicated for operation at ambient temperature and depending on the model for pressures ranging from atmospheric to UHV.
- The ECS drive series features implemented crossed roller bearings and is thus specified for high loads of up to several kilograms and guiding errors of less than 0.1mrad in pitch, yaw, and roll.
- This powerful performance is supplemented by travel ranges up to 50mm, step sizes as small as 50nm, and optional position sensor for closed loop operation with 1nm resolution. Attocube's Industrial Line positioners are available in a wide variety of designs, sizes, and travel ranges and can be stacked directly on top of each other for multi axis operation.

Specifications

| Exterior | |  |  |  |  |  |  |  |
|-------------------------------------|------------------------|---|---|---|--|---|---|---|
| Products Name | | X axis positioner | X axis positioner | X axis positioner | Goniometer | Goniometer | Rotator | Rotator |
| Part Number | | ECS3030 | ECS3080 | ECS5050 | ECGt5050 | ECGp5050 | ECR3030 | ECR5050 |
| Closed-loop travel properties | Position Resolution | 1nm | 1nm | 1nm | 0.000001° | 0.000001° | 0.00001° | — |
| | Position Repeatability | 50nm | 50nm | 50nm | ±0.00005° | ±0.00005° | ±0.0005° | — |
| | Accuracy | < 0.01% of travel range | < 0.01% of travel range | < 0.01% of travel range | ≐0.001° | ≐0.001° | ≐0.002° | — |
| Open-loop travel properties | Minimum step size | 50nm | 50nm | 50nm | 0.0001° | 0.0001° | 0.0004° | 0.0002° |
| | Fine positioning range | 1.6µm | 1.6µm | 1.6µm | 0.0014° | 0.0011° | 0.012° | 0.006° |
| Travel Range | | 20mm | 50mm | 30mm | 10° | 10° | 360° | 360° |
| Table Dimensions [mm] | | 30×30 | 30×80 | 50×50 | 50×50 | 50×50 | 30×30 | 50×50 |
| Positioning Slide | | Crossed Roller Bearing | Crossed Roller Bearing | Crossed Roller Bearing | Crossed Roller Bearing | Crossed Roller Bearing | Ball Bearing | Ball Bearing |
| Weight [kg] | Aluminum | 0.029 | 0.078 | 0.070 | 0.137 | 0.137 | 0.28 | 0.1 |
| | Stainless steel | 0.051 | 0.147 | 0.247 | 0.247 | 0.247 | 0.66 | 0.215 |
| Maximum Travel Speed | | 4.5mm/sec | 4.5mm/sec | 4.5mm/sec | ≐1°/sec | ≐1°/sec | ≐10°/sec | ≐5°/sec |
| Load Capacity (horizontal mounting) | | 9kg | 24kg | 15kg | 1kg | 1kg | 2kg | 2kg |

* Closed-loop properties are available if "encoder option" is selected. For ordering, add suffix "/NUM" after the model number.

Control of Industrial Line Positioners | ECC100

The three axes controller ECC100 is used for driving all ECS positioners either in open loop or closed loop mode, depending on the corresponding positioner model.



Software features

The ECC100 is delivered with a basic software package providing windows based software including LabView driver set and DLL. A dedicated software package /PRO providing enhanced controller functionality can be upgraded at any time. The /SYNC option offers the use of an Ethernet interface and allows for using the controller with Epics drivers.

Controller Hardware

| Part Number | ECC100 |
|-------------------|-------------------------------|
| Chassis [mm] | about (W)210 × (H)50 × (D)210 |
| Weight [kg] | 1.9 |
| Power supply | 100/115/230V, 50–60Hz |
| Power consumption | max. 100W |

Output Signals

| | |
|---------------------------------|---|
| Stepping - voltage range | 0 – 45V |
| Stepping - frequency range | 0 – 5kHz (1 axis) 0 – 2kHz (3 axes simultaneously) |
| Stepping - maximum current | >5A Peak |
| Resolution of signal generation | 680µV (16 bit) |
| Output connectors | SubD 15 connectors |

Application Systems

Optics & Optical Coatings

Holders

Bases

Manual Stages

Actuators

Motorized Stages

Light Sources

Index

Guide

Controllers/Drivers

Softwares

Stepping Motor

AC Servo Motor

Cables

Piezo

X Translation

Theta Rotation

Goniometer

Vacuum

Options

 40mm 60mm 80mm 85mm 100mm 120mm

Others